




TECHNOLOGY IN MOTION



*Baropodometria
Stabilometria
Biofeedback
Goniometria
Videografica morfologica*

www.sensormedica.com



Sensor Medica was founded in 2011 after twenty years of experience in the areas of foot pressures & gait, biomechanics, posture and CAD-CAM orthotic systems.

We study, design and develop entirely in Italy high-tech products and software for the medical-scientific and sports fields, offering a real evolution in diagnostic systems, for the optimization of human biomechanical balance.

Sensor Medica technologies are actually available in 48 countries with over **8,000** software licenses purchased and **500** CNC milling systems working.

Thanks to a commitment to ongoing research, we can produce sensor systems among the most advanced of the sector: from the most evolved detection systems, to the most complete analysis software, all completely integrated with cad-cam for orthotic production.

Sensor Medica works constantly to update technology and products to ensure them among the most advanced and innovative in the world.

freeMed Pressure sensor platforms

RunTime Pressure sensor treadmill

Moover 3D inertial motion sensor

FlexinFit In-Shoe bluetooth sensor insoles

Optoelectronic 2D-3D foot scanners, HD cameras

freeStep Biomechanics software

easyCAD Insole Orthotic modelling software

Vulcan Vx1 - twinCAM CNC milling machines for orthotic production

freeMed®

freeMed platforms are produced in aluminium for detecting plantar pressures, gait cycle with spatiotemporal measurements and balance during standing and walking.

They are lightweight and transportable, always complete of two passive plates.

High sample frequency of more than **400Hz** in real time.

All resistive sensors are **24K** gold coated to ensure an extreme reliability and repeatability and have a durability of more than one million work cycles.

All acquired information is processed in **freeStep** analysis software.

Maintenance free with no hardware re-calibration required.


3 year warranty.

freeMed® highlights

Resistive sensors, 24 K gold coated, conductive rubber foam

Sampling frequency up to 500Hz in real time; USB 2.0 interface

Bluetooth option with rechargeable batteries, sampling rate over 100 HZ



RunTime®

RunTime treadmill has **4,800** resistive sensors under the carpet for detecting plantar pressure, gait cycle with spatiotemporal measurements and balance during standing, walking and running. All sensors are **24K** gold coated to ensure an extreme reliability and repeatability.

High sample frequency of more than **200Hz** in real time.

Normally used with two synchronized cameras for studying the body posture.

All acquired information is processed in our freeStep analysis software. 3D visualization of the steps, immediate evaluation of asymmetries, curves and graphs, angle measurements on video.

3 year warranty.

RunTime® highlights

Motor 2.5 HP AC - 4 HP peak

Speed: 0.5 - 14 mph

Inclination: 15 %

Sensorized area: 21 x 59 in

Max user weight: 330 lbs

Weight: 265 lbs

Dimensions: 81 x 35 x 62 in

Dimensions folded: 55 x 35 x 64 in



Moover®

Moover is an inertial motion (IMU) sensor able to evaluate and measure your motions, accelerations and revolutions in space.




Moover® highlights

Its application in the scientific sphere allows the goniometric evaluation (range of movement) of any articulated joint. Its measurements are very important both for prevention as well as for rehabilitation.

It's also useful for studying mobility and for monitoring and documenting improvements during therapy with the ability to compare different examinations with the normal values. Moover is small, light, wireless, extremely precise and equipped with a long life battery.


freeStep software is easy to use, converts the mechanical motions into an electrical signal and always gives you normal values for comparison and generates an automatic test report.



Range of Motion



Cervical Shoulder Elbow Back Hip Knee Ankle



FlexInFit®

FlexInFit has more than 400 resistive sensors.

The system allows you to perform precise analysis inside a patient's shoe. Capable of streaming live or recording the data via SD card.

It is a versatile tool with many applications for many types of professionals: foot specialist who wants to integrate in-shoe data into his gait analysis system, physical therapist who wants to check the results of therapy, athletic trainer interested in the study and improvement of sports movements, physician interested in verifying the real pressure points inside the shoe to prevent the formation of ulcers in diabetic patients.

It can also be used to determine the efficacy of custom orthotic correction.

It is totally wireless and junction box free, to avoid any kind of interference.

FlexInFit® highlights

- Bluetooth data transmission up to 100 meters in open space
- Li-Po batteries, USB charge
- Battery life up to 4 hours
- Real time sampling frequency from 25 to 50 Hz on streaming - 250 Hz SD card
- Auto identification of foot size
- 214 resistive sensors on each insole
- Measurement scale: 0-100N/cm²
- Sensitivity: 0.1N/cm²
- 10 Bit resolution
- Thickness 0,3mm
- Material: polyester
- Very flexible, can be cut to size



Optoelectronic Systems

Podoscan 3D is a high definition laser device for acquiring the foot morphology in 3D and also showing a 2D image.

It has a precision of 1mm, allows the acquisition the foot image in neutral, semi weight bearing and weight bearing positions. It can also scan foam boxes, resin slipper socks or plaster molds.

Podoscan 2D is making a digital analysis of footprints and plantar loads. freeStep software allows to capture and store feet images and automatically create measures, comparing asymmetries and foot angles.

High-definition cameras to analyse posture and movements; video systems synchronized to our freeStep software allow a complete morphological study of the patients.



freeStep®

freeStep is an advanced software for the study of plantar pressures, motion analysis, posturography, biomechanics and humans-space relationship. It can perform many types of biomechanical evaluation including static and dynamic plantar pressures, motion analysis, EMG, morphological video analysis, digital foot scanners, posturographic studies with auto generated reports.

- Static, dynamic and balance analysis
- 2D - 3D foot scanner acquisition
- EMG detecting
- Pressure motion treadmill analysis
- One touch print and comprehensive reports
- 2D videography for asymmetry's evaluation
- Biofeedback for rehabilitation and muscles strengthening
- Statistics
- Live updates

PRESSURES analysis

Pressure sensor devices directly measure:

1. Pressures
2. Surfaces
3. Times

Pressure sensor devices indirectly measure:

1. Ground reaction force
2. Speed
3. Angles

PRESSURE analysis

What foot plantar pressure is?




$$p = \frac{F_{\perp}}{S}$$

Pressure is an intensive physics measure that is defined as the ratio between the orthogonal force module on a specific surface and on its area.




PRESSURE analysis

What foot plantar pressure is?



PRESSURE analysis

- The pressure is not directly connected with the weight of the patient
- The pressure indicates the zones that has got much more load on the tissue
- During gait the pressure (maximum) indicates the real zone where the foot goes through
- The pressure can drive the process of construction of insoles
- In an advanced point of view, the pressure can be also useful in studying posture and muscular chains function



PRESSURE SENSOR PLATFORMS

PROS

Natural walk

Self-selected speed

Barefoot / Shoes

CONS

Limited step number

Self-selected speed

Acceleration / Deceleration

PRESSURE SENSOR TREADMILL

PROS

High step number

Constant speed

Barefoot / Shoes

Running

CONS

Unknown walk way

Constant speed

Hight from the floor

Carpet vs Sensors

Not a natural Walk

IN-SHOE SENSOR INSOLES

PROS

Freedom in movement

Self-selected speed

Shoes

Shoes vs Diabetic Patients

*Feedback from corrective
elements on orthotics*

CONS

Self-selected speed

Shoes

Barefoot

STATIC analysis on platforms

1. Invite the patient on the platform,
2. Start the acquisition,
3. Make sure that his feet are aligned on same horizontal axis,
4. Position his feet in comfort position and start acquisition.
5. For balance control, repeat the static acquisition with the feet open on 30 degrees and the heels divided by 1 inch.

It offers information about the effects on foot loading of postural alignment.


STATIC analysis on platforms

- Loads and global footprints
- Maximum load pressure points
- CoP and CoF alignment
- Latero-lateral and anterior-posterior balance
- Rear-forefoot relations
- Geometric analysis
- Surface distribution



STATIC analysis on platforms


Morphological and numerical evaluation




Focus on:

- Isthmus presence
- Load points
- Adduction foot
- Width between the points
- Load on the heel and midfoot
- Aligning pressure centers
- balancing
- Rearfoot-forefoot relation

Geometric evaluation of polygon, analysis of pressure centers and loading areas



	Left foot	Right foot
Podalic angle °	14°	15°
Podalic axis °	11°	8°
Misalignment C	2.21 cm - LF A	
Misalignment L-R	0°	




	Forefoot		Middlefoot		Rearfoot	
	A	B	C	R	E	F
Surface LF cm ²	26	36	32	8	20	22
	62		40		42	
Surface RF cm ²	25	34	25	2	17	19
	59		27		36	
Load LF %	9	11	7	1	10	14
	20		8		24	
Load RF %	9	10	6	1	8	11
	19		7		19	

Pay attention to:

- CoP feet axis
- rotation of the centers of pressure
- areas with high or low pressure and areas with high or low loads

Automatic Analysis Report



Static analysis performed on the patient shows the following:

center of gravity C is located in front with respect to the polygons total support.

Centers of pressure of the legs are disaligned. Respect to the body center of gravity, the left one is in the back, the right is good.

Point of maximum pressure M is located, according to the rule in the rear right area. Load distribution between left and right leg is in the standard value (50% of the weight to the left, 50% right).

Load distribution between forefoot and rear foot, both right and left, is in the physiological values.

The surfaces of the two feet are quite similar to each other (110cm² to the left, 104cm² to the right).

Between the two rearfoot is detectable a slight difference of surface, more in the right. Between the two forefoot is detectable an excessive difference of surface, more in the left.

Comparison with normal values

Surface cm²
87

Load %
25

Weight ratio R/F
47



Surface cm²
75

Load %
23

Weight ratio R/F
50

Surface cm²
59

Load %
29

Weight ratio R/F
53



Surface cm²
51

Load %
23

Weight ratio R/F
50

146	Surface cm ²	126
54	Load %	46
521	P. Max gr/cm ²	521
218	P. Avg. gr/cm ²	215
14	Podalic angle °	15
11	Podalic axis °	8

Surface cm²
60

Load %
23

Weight ratio R/F
45 ± 3



Surface cm²
60

Load %
23

Weight ratio R/F
45 ± 3

Surface cm²
50

Load %
27

Weight ratio R/F
55 ± 3



Surface cm²
50


Load %
27

Weight ratio R/F
55 ± 3

110	Surface cm ²	110
50 ± 3	Load %	50 ± 3
-	P. Max gr/cm ²	-
-	P. Avg. gr/cm ²	-
-	Podalic angle °	-
9-15	Podalic axis °	9-15



STATIC analysis

Ideal condition




STATIC analysis

Calcaneus varus and foot supination



STATIC analysis


Calcaneus valgus and foot pronation




DYNAMIC analysis on platforms

1. Start acquiring patient data.
2. Ask the patient to walk without shoes on the platform. Make him walk at least 2-3 times, to get used to the contact with the platform.
3. Start the motion acquisition for capturing at least 8-10 footsteps.
4. You can repeat the acquisition with shoes and with shoes and custom insoles.

It offers information about the function of the muscular chains on foot loading during a dynamic movement




Biomechanical analysis, Gait line, Max line



Pay attention to:

Shape of the feet.
Areas of higher load.
Abduction of feet.
Width between the footsteps.
Movement of the CoF (white line).
Variability in footprints.

Global assessment



Pay attention to:



1. Shape of the feet
2. Movement of the CoF line (white line)
3. Areas of changing of direction
4. Movement of the red line (maximum pressure points)
5. Pressure between the 1st and 2nd toe
6. Pressure distribution between right and left

Progression of CoF (Gait Line)

The center of gravity projected on the ground (center of pressure) is considered ideal if the support follow this trend:

- a) starting centrally on the heel
- b) moving forward in the direction of V
- c) then goes on IV
- d) continues on III
- e) then on II
- f) ends on the big toe

Curves, graphs, geometric analysis, medial-lateral load distribution



Pay attention to:

1. Foot shape
2. CoF length (Gait line)
3. CoF progression should not medialize
4. Index of flex-extension



Global CoF index

Start CoF index


End CoF index

Angolo podalico


CoF progression




Normal




Flat Foot




Supinated




Pronated




Cavus - Varus




Cavus - Talus




Metatarsalgia



Body mass shifted
forward




Instability




Clubfoot

Gait cycle

Pay attention to:
1. time of double support
2. time of oscillation




Rockers



Analysis of the step rockers allows a detailed study of the four main steps of pressure distribution:

- Heel rocker
- Ankle rocker
- Forefoot rocker
- Big toe rocker

Avarage footprints




The average is a mathematic overlap of all the footprints. Are not the reality of the walking, but they offer a good point of view on the variability of the steps.

Pay attention to:

- Contours to the heel
- Contours of metatarsal area
- Contours of toes


Dynamic Analysis on platform

Comparison with normal values



Dynamic Analysis on platform

Automatic Analysis Report



Automatic pressure report
DYNAMIC ANALYSIS

The dynamic examination performed on the patient Paziente Modella the following is noted:

The total area of the left foot (98 cm²) is lower than the right (120 cm²). The length of the left footprint is 265 mm, that of the right of 270 mm.

The average distribution of the load shows appreciable differences between the left and right foot (average pressure LF = 720 gr/cm² - RF = 634 gr/cm²). The point of maximum pressure (P.Max) on the left foot is 1996 gr/cm² and on the right is 2844 gr/cm².

The load distribution between forefoot and rear foot, in the left foot, is within the physiological values (forefoot = 60% - rear foot = 40%), while the right shows a rear-forefoot distribution not conform to normal (forefoot = 70% - rear foot = 30%). The lateral distribution of the load, on the left foot, is 43% on lateral and 57% on medial, on the right side is 47% on lateral and 53% on medial.

The stay on the ground of the left foot was measured in 956ms for a total of 154 frames, for the right 1022ms in 179 frames. During the frame 116 of the left foot there is the point of maximum load, whereas the right took place during the frame 137.

Pressure detection during dynamic shows the following:

feet surfaces are similar to each other (SX = 109 cm² - DX = 106 cm²).

Length of the left foot is 230 mm, right is 221 mm.

Load average distribution shows differences between the left foot and the right (Left = 974 g/cm² - Right = 894 g/cm²).

Load point (P Max) in the left foot is 2684 g/cm² and in the right is 2448 g/cm².

Left load distribution between forefoot and rearfoot is in the physiological values (forefoot = 57% - rearfoot = 43%), while the right shows a nonconformity distribution with the standard values (forefoot = 55% - rearfoot = 45%).

Latero-lateral load distribution for the left foot is 52% in lateral and 48% in medial, for right foot is 38% in lateral and 62% in medial.

Speed of the left foot is 723ms for a total of 153 frames, right is 664ms in 155 frames.

Step 104 of left foot registers the maximum load, while step 45 of the right foot.

BALANCE analysis on platform

1. Place the patient on the platform
2. Put patient data inside the software
3. Make sure that the feet are aligned on the horizontal axis
4. Open the patient feet of 30 degrees and separated the heels of 2 cm.
5. Select the acquisition protocol and begin the analysis

Protocols of Posturographic analysis

Selected test

Sway Test

Sway Test

Romberg Test

Cervical Test

Stomatognathic Test

Postural test

Protocollo Monopodalico Funzionale

Protocollo ATM

Selected test

Sway Test

(Freq. 50Hz)

Test

Duration

Bipodalic OE

51,2



Bipodalic CE

51,2



Monopodalic LF OE

10,0



Monopodalic RF OE

10,0



Monopodalic LF CE

10,0



Monopodalic RF CE

10,0




Stability test: measure COP movements


The COP trace is shown in 2D (left-right on X-axis, back-forward on Y axis) for explaining the patient instability during the test.

The ellipse automatically include the 90% of the total trace.

Ellipse surface gives information on the efficiency and the balancing of the tonic-postural system.



Stability test: measure COP movements central, left and right



The COP tracking, in foot pressure platforms, is the result of right lateral and left lateral center of gravity swing.

With this evaluation, we can study not only the global postural attitude but also the dynamic behaving of the two legs.

STABILITY test

Stabilometric analysis measures the COP path in function of the time on X-Y coordinates.

Latero-lateral and anterior-posterior movements.

All measurements are mathematical calculations and geometric derivatives



RunTime

Pressure Sensor Treadmill

Acquisition techniques and main indications

Dynamic analysis on RunTime treadmill

1. Place the patient on the treadmill
2. Put the patient data in the software
3. Start slowly the treadmill
4. Increase the speed until you reach the desired speed
5. Let the patient run or walk for at least 2 minutes
6. Start the acquisition, normally for minimum 30 sec


Total Test duration s	10,92		
Number of analyzed footprints	15		
Average speed Km / h	2,16		
Cadence Steps per min	82		
	LF	RF	Simmety Index
Length Half/Step cm	44,94	47,25	95%
Stance time ms	927	1009	92%
Dual load/flight Time ms	244	269	90%
Rocking Time ms	527	502	95%
Rockers ICP MS	134	151	89%
Rockers FFCP ms	520	578	90%
Rockers FFP ms	262	277	95%
Surface cm ²	192	203	95%
P.max gr / cm ²	838	755	90%
Average P. gr / cm ²	806	734	91%
Fore-foot Load %	82	81	99%
Rear-foot Load %	18	19	97%
Medial Load %	41	47	89%
Lateral Load%	59	53	91%
Podalic Axis °	6,38	3,86	61%

Avarage global value

Automatic comparison between left and right steps are displayed with percentages index.

Green indicators display a good balance inside the range of 10%

From orange to red there is an incremental asymmetry.



Biomechanical Analysis


- Time of left and right foot support, time of double support.
- Butterfly chart, for the study of the centers of dynamic pressure.
- 3D animation in slow motion allows to analyze frame-by-frame the distribution phases synchronized with cameras.



Indicators

Trend of the main values of feet support

Comparison between right and left



Gait-line of each single step

1. Medial-lateral loads
2. Length of footprint and of the gait line
3. Feet angles and axis


FlexInFit

Bluetooth Sensor insoles


Acquisition techniques and main indications

Dynamic Analysis with FlexInFit

1. Select the proper insoles size
2. Insert the insoles in the shoes and connect the electronic transmitters
3. Switch on the 2 device, wait the activation of Bluetooth connection
4. Let the patient walk naturally
5. Start the acquisition, normally for minimum 30 sec



Gait-line and Butterfly




Videography

Acquisition techniques and main indications

2D Videography

Videography allows measuring the morphological posture of the patient for the identification of asymmetries in static and in dynamic.

An automatic protocol guide the positioning of markers for the calculation of lengths and angles of the body segments.




2D VIDEOGRAPHY: morphological asymmetries


2D videography allows the complete morphological measurement of the subject for the identification of body asymmetries.

The analysis is facilitated by a guided marker protocols.

Also free measurements are allowed.



2D VIDEOGRAPHY: morphological asymmetries



Point

Apex of the patient
C7 vertebra
Right shoulder
Left shoulder
Right scapula
Left scapula
Right sips
Left sips
Right wrist
Left wrist
Trochanter right
Trochanter left
Right popliteal
Left popliteal
Right Achilles tendon
Left Achilles tendon
Right heel
Left heel


Results

Shoulder tilt	1° LF^
Ground distance left shoulder mm	1369
Ground distance right shoulder mm	1373
Scapulae tilt	4° LF^
Sips tilt	3° LF^
Angle varus/valgus knee rg	164°
Angle varus/valgus rearfoot rf	176°
Angle varus/valgus knee lf	165°
Angle varus/valgus rearfoot lf	177°

2D Videography

Start the motion capture during the walking.

Image capture in static is instantaneous, in dynamic is automatically stopped at the end of the movement.




2D PodoScan

Acquisition technology and main indications

1. Place the patient on the scanner, with both feet on the glass
2. Start the capture interface of the program
3. Wait until the scanning is complete (2-3 sec)
4. Get out patient from the scanner






Automatic calculations

Foot length
 Foot size
 Forefoot width
 Arch length
 Heel width
 Heel length
 Angle of the feet
 Angle of the metatarsals


Other automatic measurements



A protocol wizard facilitates the measurement of the foot. Measurement free hand tools ensure full independence

Different images

1. Natural image
2. Color Image to highlight contact areas
3. Gray image that recalls the podographic paper




3D Laser PodoScan

Acquisition technology and main indications

Acquisition

1. Place the subject on the scanner with the barefoot in half load
2. Start the acquisition in **freeStep**
3. Wait until **laser** scanning is completed (2 seconds)



Collecting data for orthotic design

Allowing the height measurement **3D PodoScan** ensures a perfect, three-dimensional reconstruction of the foot arch, wrapping the heel and the area of the forefoot.