



Software Installation Instructions

- Download and install the demo version of freeStep from <http://win.sensormedica.com/update/freestep/setupfreestep.exe>
- Create an appointment for activation by contacting your representative or visiting support.sensormedica.com
- In a web browser, navigate to join.zoho.com and enter the session ID provided.
- Your assigned technician will activate the software and train you on use.
- Please note that software user manual is integrated into the software by pressing F1 in the software or via the drop down menu.

Computer Minimum Requirements:

- Windows 10
- i5 processor or higher
- 16 Gig Ram or higher
- 256-Bit Graphics Card
- Solid State Hard Disk (Highly Recommended)

Sensor Medica USA LLC
919 South Utah Avenue
Idaho Falls, ID 83402 USA
www.sensormedicausa.com
usa@sensormedica.com



Technology in Motion